

2009

■ GEI SV
PCR SV
CLEANUP SV
COMBO GP

expin[™]

Expin[™] Kits Protocol
Handbook

For purification of fragment DNA


GeneAll

Customer & Technical Support

Do not hesitate to ask us any question.

We thank you for any comment or advice.

Contact us at

www.geneall.com

Tel : 82-2-407-0096

Fax : 82-2-407-0779

E-mail(Order/Sales) : sales@geneall.com

E-mail(Tech. Info.) : bnp@geneall.com

Visit GeneAll® Community

www.geneall.com

www.geneall.co.kr

Trademarks

Exprep™, Exfection™, Expin™, Exgene™, EzSep™, EzClear™, RiboEx™
Ribo_clear™ and GenEx™ are trademarks of GeneAll Biotechnology co., Ltd.

© 2009 GeneAll® Biotechnology, all right reserved.

This protocol handbook is included in :

GeneAll® Gel SV (102-150, 102-102)

GeneAll® PCR SV (103-150, 103-102)

GeneAll® CleanUp SV (113-150, 113-102)

GeneAll® Combo GP (112-150, 112-102)

Visit www.geneall.com or www.geneall.co.kr for FAQ, QnA and more information.

→ INDEX

↳ GeneAll® Expin™ Protocols

Kit Contents	4
Storage Conditions	5
Quality Control	5
Safety Information	5
Introduction	6
Product Specifications	7
Applicable range in length of DNA	7
General considerations	8
Binding	8
Washing	9
Elution	9
Correlation between the elution volume and the recovery	10
Electrophoresis analysis	11
Protocols for..	
Gel SV Protocol	12
PCR SV Protocol	18
CleanUp SV Protocol	22
Troubleshooting Guide	26
Ordering Information	28
Expin™ Gel SV Brief protocol	30
Expin™ PCR SV Brief protocol	30
Expin™ CleanUp SV Brief protocol	31

Kit Contents

Cat. No.	Gel SV		PCR SV	
	102-150	102-102	103-150	101-102
Size	mini	mini	mini	mini
No. of preparation	50	200	50	200
GeneAll® SV Column type D	50	200	50	200
Collection Tube	50	200	50	200
Buffer GB	60 ml	120 ml x 2	-	-
Buffer PB	-	-	30 ml	120 ml
Buffer NW	60 ml	250 ml	60 ml	250 ml
Buffer EB*	15 ml	30 ml	15 ml	30 ml
Protocol Handbook	1	1	1	1

Cat. No.	CleanUp SV		Combo GP	
	113-150	113-102	112-150	112-102
Size	mini	mini	mini	mini
No. of preparation	50	200	50	200
GeneAll® SV Column type D	50	200	50	200
Collection Tube	50	200	50	200
Buffer GB	-	-	60 ml	120 ml x 2
Buffer PB	-	-	30 ml	120 ml
Buffer NR	30 ml	120 ml	-	-
Buffer NW	60 ml	250 ml	60 ml	250 ml
Buffer EB*	15 ml	30 ml	15 ml	30 ml
Protocol Handbook	1	1	1	1

* 10mM TrisCl, pH 8.5

Storage Conditions

All components of GeneAll® Expin™ Gel, PCR and CleanUp SV Kit should be stored at room temperature (15~25°C) and are stable for one year. Protect from exposure to direct sunlight.

During shipment or storage under cold ambient condition, a precipitate may be formed in buffer GB, PB or NR. In such a case, heat the bottle at 50°C to dissolve completely. Using precipitated buffers will lead to poor DNA recovery. GeneAll® Expin™ SV kit series are guaranteed until the expiration date printed on the product label.

Quality Control

All components in GeneAll® Expin™ Gel, PCR and CleanUp SV Kit are manufactured in strictly clean condition, and its degree of cleanness is monitored periodically. Restriction enzyme assay, PCR amplification assay, spectrophotometric assay and automatic sequencing as the validation of quality are carried out from lot to lot thoroughly, and only the qualified is approved to be delivered.

Safety Information

The buffers included in GeneAll® Expin™ Gel, PCR and CleanUp SV Kit contain irritant which is harmful when in contact with skin or eyes, or when inhaled or swallowed. Care should be taken during handling. Always wear gloves and eye protector, and follow standard safety precautions.

Buffer GB, PB and NR contain chaotropes. It can form highly reactive compounds when combined with bleach. Do NOT add bleach or acidic solutions directly to the sample-preparation waste.

Introduction

GeneAll® Expin™ SV series provide reliable and fast methods for the purification of fragment DNA from agarose gel and PCR/enzyme reaction mixtures. GeneAll® Expin™ SV series consist of;

GeneAll® Expin™ Gel SV (102-150, 102-102)

Gel SV kit is designed for a fast and efficient isolation of 80 bp to 10 kb of DNA fragments from standard or low-melting agarose gel in TAE or TBE buffer system. Purified DNA can be directly used in ligation, labelling, sequencing and many other downstream application without further manipulation.

GeneAll® Expin™ PCR SV (103-150, 103-102)

PCR SV kit provides a simple and rapid method to purify PCR products or other enzymatic reactions in just 6 minutes. Up to 10 ug of pure DNA (100 bp~ 10 kb) can be obtained, and this purified DNA can be directly used in cloning, sequencing and many other application. PCR SV procedures remove the DNA fragment smaller than 100 bp, resulting in removal of primers and primer dimers in PCR products.

GeneAll® Expin™ CleanUp SV (113-150, 113-102)

CleanUp SV kit provides a simple and rapid method to purify 40 bp ~ 10 kb DNA from enzymatic reactions in just 6 minutes. Up to 10 ug of pure DNA which is at least 40 bp but less than 10 kb in length can be obtained using this kit and the purified DNA can be directly used for sequencing, cloning and other routine applications without further manipulation.


GeneAll® Expin™ Combo GP (112-150, 112-102)

The combination set of GeneAll® Gel SV and PCR SV.

Product Specifications

GeneAll® Expin™	Gel SV	PCR SV	CleanUp SV
Standard sample size	200 mg of gel (400 mg max.)	100 ul of PCR reactions	50 ul of enzyme mixtures
Recovered DNA size	80 bp~10 kb	100 bp~10 kb	40 bp~10 kb
Typical yields	~ 85 %	~ 95 %	~ 95 %
Maximum binding capacity	10 ug	10 ug	10 ug
Preparation time	~ 15 mins	~ 6 mins	~ 6 mins
Minimum elution volume	30 ul	30 ul	30 ul

Applicable range in length of DNA



GeneAll® Expin™ SV series consists of Gel, PCR and CleanUp SV Kit. Each kit is optimized for efficient recovery of DNA and removal of contaminants in each specific application.

General considerations

GeneAll® Expin™ SV series takes advantage of silica membrane and spin/vacuum column technology to recover DNA fragments. Under high salt conditions DNA binds to silica membrane and impurities pass through membrane into a collection tube. The membranes are washed with an ethanol-containing buffer to remove any traces of proteins, salts, remnant of agarose and other enzymatic reaction components. Finally pure DNA is released into a clean collection tube with water or low ionic strength buffer.

Binding

The basic principle which DNA binds to silica membrane is identical in all GeneAll® Expin™ SV series. Binding buffers (Buffer GB, PB and NR) make the optimal binding condition in each specific applications. Buffer GB in Gel SV kit is composed to dissolve standard agarose gel as well as low melt agarose gel in addition to adjusting a binding condition. Usually low melt agarose gel results a better recovery yield. DNA binds to silica membrane at lower pH than pH 7.5. The components and pH of starting sample can alter the pH of the mixture with binding buffer, especially in extraction of DNA from agarose gel. Buffer GB contains pH indicator in order to check this alteration of binding condition. If the color (yellow) of binding mixture turns to brown or purple after addition of buffer GB, it means that the pH of binding mixtures is higher than the optimal, and it can be easily adjusted with small volume of sodium acetate before proceeding with the protocol. The indicator dye is completely removed during subsequent washing steps and does not interfere the downstream applications.

Washing


Any unwanted oligos and impurities, such as salts, proteins, nucleotides, agarose, dyes and detergents will not be bound but be passed through the silica membrane. A minute impurities, such as salts, are quantitatively washed away with buffer NW which contains ethanol. The quality of DNA may be slightly increased with the repeat of washing. Any residual ethanol should be removed completely with an additional centrifugation because the residual ethanol in eluate may interfere some subsequent applications.

Elution

DNA is released under the condition of low salts and neutral or weakly alkaline pH ($7 < \text{pH} < 9.5$). Although buffer EB (10mM TrisCl, pH 8.5), TE, or distilled water can be used for elution, it should be considered that EDTA in buffer TE may interfere the subsequent reactions and low pH (< 7.0) of distilled water can reduce DNA recovery. Because water does not have any buffering agents the eluate by water should be stored under -20°C not to degrade.

The minimum elution volume is 30 μl and lower volume will decrease the yield significantly. It is important for optimal elution to apply the elution buffer to the center of the membrane, because the membrane should be covered completely by eluent for an optimal recovery. Up to 200 μl of elution buffer can be applied and it results in low concentration of DNA. Higher concentrated DNA will be obtained with lower elution volume, and maximum yield can be obtained by larger elution volume. The yield with large fragments (> 5 kb) can be increased slightly by using pre-warmed (70°C) elution buffer. Incubation for 1 minute after addition of elution buffer may increase the efficiency of elution.

Correlation between the elution volume and the recovery




Correlation between the elution volume and the recover rates. Upper panel; GeneAll[®] PCR SV Kit. Lower panel; GeneAll[®] Gel SV Kit. 3 ug of 494 bp PCR products were purified and eluted with the indicated volumes of buffer EB. Elution volume lower than 30 ul causes significant loss of DNA. 1/10 volume of eluate purified with Gel SV Kit was resolved on 1 % agarose gel; lower right. (1/20 for PCR SV Kit, upper right)


Recovery rates (%)			
DNA size (bp)	Gel SV	PCR SV	CleanUp SV
60	39	0	63
120	71	78	80
200	76	83	84
800	84	94	94
1800	82	91	93
4300	78	85	88
8700	73	76	79

Average recovery rates of GeneAll® Expin™ SV Kit with various sizes of DNA. 3 ug of starting sample was purified and eluted with 50 ul of buffer EB. Optional steps were not performed and SV columns were incubated for 1 minutes after addition of buffer EB.

Electrophoresis analysis



Quantities of extracted 4.5 kb DNA fragment corresponding to 1/5 of the DNA obtained by purification from 0.5 ug starting DNA with a recovery of 90 %. Sample were run on 1 % TAE agarose gel. M; Lambda-BstP1 marker Lane 1; Total amount before extraction (0.5 ug) Lane 2; 1/5 amount after extraction [90 ng compared to known amount (10-150 ng) DNA] * Total obtained amount of DNA = 90 x 5 = 450 ng approximately (90 %)



PCR products of several sizes (70, 176, 757 and 1487 bp from left to right) were purified using GeneAll® CleanUp SV Kit. The band intensity of 1/5 amount of 20 ul PCR products (lane 1,5,9,13) was almost identical with that of 1/5 amount of 50 ul eluate. (2~4, 6~8, 10~12 and 14~16, triplicate) M; 1 kb ladder marker

GeneAll® Expin™ Gel SV Protocol

Before experiment

- * All centrifugation should be carried out at 10,000 xg above (> 12,000 rpm) at room temperature in a microcentrifuge.
- * Prepare water bath or heating block to 50 °C.
- * All solutions should be equilibrated at room temperature before procedures.
- * For large fragments (>5 kb), pre-warm buffer EB to 70 °C
- * If Buffer GB has precipitated, dissolve it completely at 50 °C before starting

Spin/Vacuum Methods for Gel Extraction

When using this kit, one of the two methods can be chosen to purify DNA. DNA fragment can be purified using centrifugation to pull the mixture through the SV column. Alternatively, vacuum can be used to force the mixture out through the column.



Centrifugation Protocol

1. Excise the DNA band of interest using a ethanol-cleaned razor blade or scalpel on a transilluminator.

Minimize gel volume by cutting the gel slice as small as possible.

Use of long wave length transilluminator and short handling time will lead to better quality of DNA. It can be critical in some experiments, such as ligation.

2. Weigh the gel slice in a microcentrifuge tube. Add 3 volumes (ul) of Buffer GB to 1 volume (mg) of gel.

For example, add 300 ul of buffer GB to each 100 mg of agarose gel slice. For > 1.5 % agarose gel, add 5 volumes of buffer GB.

3. Incubate at 50°C until the agarose gel is completely melted (5 ~ 10 min).

To help the efficient dissolving of gel, vortex the tube every 2-3 min during the incubation.

4. After the slice has dissolved completely, check that the color of the mixture is yellow (similar to Buffer GB).

If the color of the mixture becomes brown or purple, add 10 ul of 3M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. This pH adjustment is not needed for the change of color due to the ingredient of loading dye (eg. bromophenol blue, xylene cyanol)

5. (Optional :) Add 1 gel volume of isopropanol to the sample and vortex to mix.

For 100 mg of gel volume, add 100 ul of isopropanol. Do NOT centrifuge at this step.

This step is required to increase the recover yields of DNA fragments <200 bp or >5 kb. For the DNA fragments between 200 bp and 5 kb, it has little effect on the recovery.

6. Transfer the mixture to a SV column. Centrifuge for 1 min. Discard the pass-through and re-insert the SV column into the Collection Tube.

If the mixture volume is larger than 700 ul, apply the mixture twice; apply 700 ul of the mixture, spin down, discard the pass-through, re-insert empty collection tube, and repeat the step again until all of the mixture has been applied to the SV column.

7. (Optional :) Apply 500 ul of Buffer GB to the SV column. Centrifuge for 30 sec. Discard the pass-through and reinsert the SV column into the Collection Tube.

This step is for further complete removal of any traces of agarose and required only for direct use of purified DNA for very sensitive applications, such as in vitro transcription. Usually this step can be skipped for automatic sequencing or ligation

8. Add 700 ul of Buffer NW to the SV column. Centrifuge for 30 sec. Discard the pass-through. Reinsert the SV column into the Collection Tube.

If the purified DNA will be used for salt sensitive applications, let the SV column stand for 5 min after addition of buffer NW, making some amount of wash buffer flow through the column by gravity before centrifugation.

9. Centrifuge for additional 1 min to remove residual wash buffer. Transfer the SV column to a new 1.5 ml tube.

If residual ethanol remains in the SV column, centrifuge again for additional 1 min at full speed before transferring to a new 1.5 ml tube. Residual ethanol from buffer NW can inhibit subsequent enzymatic reaction.

10. Apply 50 ul of Buffer EB or ddH₂O to the center of the membrane in the SV column, let stand for 1 min, and centrifuge for 1 min.

Ensure that the buffer EB or distilled water is dispensed directly onto the center of SV column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30 ul of elution buffer, but the volume lower than 30 ul will decrease the yield significantly. Up to 200 ul of elution buffer can be applied to SV column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in buffer EB (10mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

B Vacuum Protocol

The vacuum pressure should be in the range of the list below. Lower vacuum pressure may reduce DNA yield and purity.

Most commercial vacuum manifold with luer connectors can be used with this protocol.

15 ~ 18 in Hg	380 ~ 460 mbar
285 ~ 345 mm Hg	5.5 ~ 6.5 Psi

1. Excise the DNA band of interest using a ethanol-cleaned razor blade or scalpel on a transilluminator.

Minimize gel volume by cutting the gel slice as small as possible.

Use of long wave length transilluminator and short handling time will lead to better quality of DNA. It can be critical in some experiments, such as ligation.

2. Weigh the gel slice in a microcentrifuge tube. Add 3 volumes (ul) of Buffer GB to 1 volume (mg) of gel.

For example, add 300 ul of buffer GB to each 100 mg of agarose gel slice. For > 1.5 % agarose gel, add 5 volumes of buffer GB.

3. Incubate at 50 °C until the agarose gel is completely melted (5 - 10 min).

To help the efficient dissolving of gel, vortex the tube every 2-3 min during the incubation.

4. After the slice has dissolved completely, check that the color of the mixture is yellow (similar to Buffer GB).

If the color of the mixture becomes brown or purple, add 10 ul of 3M sodium acetate, pH 5.0, and mix. The color of the mixture will turn to yellow. This pH adjustment is not needed for the change of color due to the ingredient of loading dye (eg. bromophenol blue, xylene cyanol)

5. (Optional :) Add 1 gel volume of isopropanol to the sample and vortex to mix.

For 100 mg of gel volume, add 100 ul of isopropanol. Do NOT centrifuge at this step.

This step is required to increase the recover yields of DNA fragments <200 bp or >5 kb. For the DNA fragments between 200 bp and 5 kb, it has little effect on the recovery.

6. Attach the SV column to a port of the vacuum manifold tightly.

Most commercial vacuum manifold with luer connectors can be adopted to this protocol

7. Transfer the mixture to a SV column, by pipetting. Switch on vacuum source to draw the solution through the SV column. When all liquid has been pulled through the SV column, release the vacuum. Repeat this step until all mixture have applied to the SV column.

8. (Optional :) Apply 500 ul of Buffer GB to the SV column and switch on vacuum source. When all liquid has been pulled through the SV column, release the vacuum.

This step is for further complete removal of any traces of agarose and required only for direct use of purified DNA for very sensitive applications, such as in vitro transcription. Usually this step can be skipped for automatic sequencing.

9. Apply 700 ul of Buffer NW and switch on vacuum source. When all liquid has been pulled through the SV column, release the vacuum. Transfer the SV column to a empty collection tube (provided).

If the purified DNA will be used for salt sensitive applications, let the SV column stand for 5 minutes after addition of buffer NW before applying vacuum.

10. Centrifuge for additional 1 min to remove residual wash buffer. Transfer the SV column to a new 1.5 ml tube.

If the SV column has buffer NW associated with it, centrifuge again for additional 1 min at full speed before transferring to a new 1.5 ml tube. Residual ethanol from buffer NW can inhibit subsequent enzymatic reaction.

11. Apply 50 ul of Buffer EB or ddH₂O to the center of the membrane in the SV column, let stand for 1 min, and centrifuge for 1 min.

Ensure that the buffer EB or distilled water is dispensed directly onto the center of SV column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30 ul of elution buffer, but the volume lower than 30 ul will decrease the yield significantly. Up to 200 ul of elution buffer can be applied to SV column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in buffer EB (10mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

GeneAll® Expin™ PCR SV Protocol

Before experiment

- * All centrifugation should be carried out at 10,000 xg above (> 12,000 rpm) at room temperature in a microcentrifuge.
- * All solutions should be equilibrated at room temperature before procedures.
- * For large fragments (>5 kb), pre-warm buffer EB to 70 °C

Spin/Vacuum Methods for PCR Purification

When using this kit, one of the two methods can be chosen to purify DNA. DNA fragment can be purified using centrifugation to pull the mixture through the SV column. Alternatively, vacuum can be used to force the mixture through the column.



Centrifugation Protocol

- 1. Add 5 volumes of Buffer PB to 1 volume of the sample and mix. Transfer the mixture to a SV column.**
For 100 ul reaction, add 500 ul of buffer PB. It is not necessary to remove mineral oil.
- 2. Centrifuge for 30 sec. Discard the pass-through and reinsert the SV column back into the same tube.**
- 3. Apply 700 ul of Buffer NW. Centrifuge for 30 sec. Discard the pass-through and reinsert the SV column back into the collection tube.**

4. Centrifuge for an additional 1 min to remove residual wash buffer. Transfer the SV column to a new 1.5 ml tube.

If the SV column has buffer NW associated with it, centrifuge again for additional 1 min at full speed before transferring to the new 1.5 ml tube. Residual ethanol from buffer NW can inhibit subsequent enzymatic reaction.

5. Apply 50 ul of Buffer EB or ddH₂O to the center of the membrane in the SV column, let stand for 1 min, and centrifuge for 1 min.

Ensure that the buffer EB or distilled water is dispensed directly onto the center of SV column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30 ul of elution buffer, but the volume lower than 30 ul will decrease the yield significantly. Up to 200 ul of elution buffer can be applied to SV column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in buffer EB (10mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

B Vacuum Protocol

The vacuum pressure should be in the range of the list below. Lower vacuum pressure may reduce DNA yield and purity.

Most commercial vacuum manifold with luer connectors can be used with this protocol.

15 ~ 18 in Hg	380 ~ 460 mbar
285 ~ 345 mm Hg	5.5 ~ 6.5 Psi

- 1. Attach the SV column to a port of the vacuum manifold tightly.**
- 2. Add 5 volumes of Buffer PB to 1 volume of the sample and mix. Transfer the mixture to a SV column by pipetting.**
For 100 ul reaction, add 500 ul of buffer PB.
It is not necessary to remove mineral oil.
- 3. Switch on vacuum source to draw the solution through the SV column. When all liquid has been pulled through the SV column, release the vacuum.**
- 4. Apply 700 ul of Buffer NW and switch on vacuum source. When all liquid has been pulled through the SV column, release the vacuum. Transfer the SV column to a collection tube (provided).**
- 5. Centrifuge for 1 min to remove residual wash buffer. Transfer the SV column to a new 1.5 ml tube.**

If the SV column has buffer NW associated with it, centrifuge again for additional 1 min at full speed before transferring to the new 1.5 ml tube. Residual ethanol from buffer NW can inhibit subsequent enzymatic reaction.

6. Apply 50 ul of Buffer EB or ddH₂O to the center of the membrane in the SV column, let stand for 1 min, and centrifuge for 1 min.

Ensure that the buffer EB or distilled water is dispensed directly onto the center of SV column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30 ul of elution buffer, but the volume lower than 30 ul will decrease the yield significantly. Up to 200 ul of elution buffer can be applied to SV column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in buffer EB (10mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

GeneAll® Expin™ CleanUp SV Protocol

Before experiment

- * All centrifugation should be carried out at 10,000 xg above (> 12,000 rpm) at room temperature in a microcentrifuge.
- * All solutions should be equilibrated at room temperature before procedures.
- * For large fragments (>5 kb), pre-warm buffer EB to 70°C

Spin/Vacuum Methods for DNA Clean-up

When using this kit, one of the two methods can be chosen to purify DNA. DNA fragment can be purified using centrifugation to pull the mixture through the SV column. Alternatively, vacuum can be used to force the mixture through the column.

A Centrifugation Protocol

- 1. Add 10 volumes of Buffer NR to 1 volume of the sample and mix. Transfer the mixture to a SV column.**
For 50 ul reaction, add 500 ul of buffer NR.
If the length of DNA is longer than 100 bp, add 5 volumes of buffer NR.
- 2. Centrifuge for 30 sec. Discard the pass-through and reinsert the SV column back into the collection tube.**
- 3. Apply 700 ul of Buffer NW. Centrifuge for 30 sec. Discard the pass-through and reinsert the SV column back into the collection tube.**

4. Centrifuge for an additional 1 min to remove residual wash buffer. Transfer the SV Column to a new 1.5 ml tube.

If the SV column has buffer NW associated with it, centrifuge again for additional 1 min at full speed before transferring to the new 1.5 ml tube. Residual ethanol from buffer NW can inhibit subsequent enzymatic reaction.

5. Apply 50 ul of Buffer EB or ddH₂O to the center of the membrane in the SV column, let stand for 1 min, and centrifuge for 1 min.

Ensure that the buffer EB or distilled water is dispensed directly onto the center of SV column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30 ul of elution buffer, but the volume lower than 30 ul will decrease the yield significantly. Up to 200 ul of elution buffer can be applied to SV column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in buffer EB (10mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

B Vacuum Protocol

The vacuum pressure should be in the range of the list below. Lower vacuum pressure may reduce DNA yield and purity.

Most commercial vacuum manifold with luer connectors can be used with this protocol.

15 ~ 18 in Hg	380 ~ 460 mbar
285 ~ 345 mm Hg	5.5 ~ 6.5 Psi

- 1. Attach the SV column to a port of the vacuum manifold tightly.**
- 2. Add 10 volumes of Buffer NR to 1 volume of the sample and mix. Transfer the mixture to a SV column by pipetting.**
For 50 ul reaction, add 500 ul of buffer NR.
If the length of DNA is longer than 100 bp, add 5 volumes of buffer NR.
- 3. Switch on vacuum source to draw the solution through the SV column. When all liquid has been pulled through the SV column, release the vacuum.**
- 4. Apply 700 ul of Buffer NW and switch on vacuum source. When all liquid has been pulled through the spin column, release the vacuum. Transfer the SV column to a collection tube (provided).**
- 5. Centrifuge for an additional 1 min to remove residual wash buffer. Transfer the SV column to a new 1.5 ml tube.**
If the SV column has buffer NW associated with it, centrifuge again for additional 1 min at full speed before transferring to the new 1.5 ml tube. Residual ethanol from buffer NW can inhibit subsequent enzymatic reaction.

6. Apply 50 ul of Buffer EB or ddH₂O to the center of the membrane in the SV column, let stand for 1 min, and centrifuge for 1 min.

Ensure that the buffer EB or distilled water is dispensed directly onto the center of SV column membrane for optimal elution of DNA.

To obtain more concentrated DNA solution, apply 30 ul of elution buffer, but the volume lower than 30 ul will decrease the yield significantly. Up to 200 ul of elution buffer can be applied to SV column, and it will reduce the concentration of DNA.

For larger fragment (>5 kb), use pre-warmed (70°C) elution buffer for best efficiency.

For long-term storage, eluting in buffer EB (10mM TrisCl, pH 8.5) or TE, pH 8.0 and storing under -20°C is recommended. When using water for elution, check the pH of water (>pH 7.0) before elution.

Troubleshooting Guide

Facts	Possible Causes	Suggestions
Low or no recovery	[Gel SV] Incompletely solubilized gel	The sliced agarose gel should be completely dissolved without any particles. To assist the complete solubilization, mix the tube by vortexing every 2-3 minutes during incubation, or increase the incubation time. To use low melt agarose usually results better recovery.
	[Gel SV] Too high pH of binding mixture	At high pH, the binding of DNA to silica membrane will be significantly reduced. The dye included in buffer GB indicates the pH of mixture as color change from yellow at optimal pH to brown or purple at abnormally higher pH. If the color of mixture has turned to brown or purple, add 10 ul of 3M sodium acetate, pH 5.0 to the sample and mix. The color of mixture will turn to yellow indicating the correct pH for DNA binding.
	Improper elution buffer	As user's requirement, elution buffer other than buffer EB can be used. The condition of optimal elution is low salt concentration with alkaline pH ($7 < \text{pH} < 9$). When water or other buffer was used as eluent, ensure that conditions.
	Elution buffer incorrectly dispensed	Ensure that elution buffer dispensed to the center of membrane. Incorrectly dispensed elution buffer causes inappropriate contact with membrane, followed by poor DNA recovery.
Ligation failure	[Gel SV] Too long or strong exposure to UV on transilluminator	UV destroys the DNA ends. Use UV of long wave length and make the handling time as short as possible when excising the gel slice.

Facts	Possible Causes	Suggestions
Column clogging	[Gel SV] Incompletely solubilized gel	See the section 'Incompletely solubilized gel' in the Facts "Low or no recovery"
	[Gel SV] > 1.5% agarose gel is used	For > 1.5 % agarose gel, 5 volumes of buffer GB to 1 volume of gel slice should be added. For 100mg of agarose gel, add 500 ul of buffer GB. If the column is clogged, transfer the mixture from the SV column to a 1.5 ml tube, add 1 volume of buffer GB to mixture volume. Incubate for 5 minutes at 50°C, proceed again to binding steps.
Enzymatic reaction is not performed well with the purified DNA	Residual ethanol from Buffer NW remains in eluate	It is essential to remove any residual ethanol included in buffer NW from column membrane. Centrifuge again for complete removal of ethanol.
	Too high salt concentration in eluate	Incubate for 5 minutes after addition of buffer NW at washing steps.
	Eluate contains denatured ssDNA	For reannealing of ssDNA to dsDNA, incubate ssDNA at 95°C for 2 min, and then allow to cool slowly to room temperature.
DNA floats out of well while loading of agarose gel	Residual ethanol from Buffer NW remains in eluate	It is essential to remove any residual ethanol included in buffer NW from SV column membrane. Centrifuge again for complete removal of ethanol.
Nonspecific band appears after purification	DNA denatured	Renature the DNA by warming up to 95°C for a minute and let cool slowly to room temperature.

Ordering Information

※ For more information about ordering, visit www.geneall.com

Products	Cat. No.	Type	Size
GeneAll® Exprep™ for preparation of plasmid DNA			
Plasmid SV min	101-150	spin/ vacuum	50
	101-102		200
	101-111		1,000
Plasmid SV Midi*	101-226	spin/ vacuum	26
	101-250		50
	101-201		100
Plasmid Quick	101-050	mini/ spin	50
	101-002		200
	101-011		1,000
GeneAll® Exfection™ for preparation of highly pure plasmid DNA			
Plasmid LE mini	111-150	spin/ vacuum	50
	111-102		200
Plasmid LE Midi*	111-226	spin/ vacuum	26
	111-201		100
Plasmid EF Midi*	121-220	spin	20
GeneAll® Expin™ for purification of fragment DNA			
Gel SV	102-150	mini/ spin/ vacuum	50
	102-102		200
PCR SV	103-150	mini/ spin/ vacuum	50
	103-102		200
CleanUp SV	113-150	mini/ spin/ vacuum	50
	113-102		200
Combo GP	112-150	mini/ spin/ vacuum	50
	112-102		200
GeneAll® Exgene™ for isolation of total DNA			
Tissue SV mini (plus!)*	104(9)-101	spin/ vacuum	100
	104(9)-152		250
Tissue SV Midi (plus!)**	104(9)-226	spin/ vacuum	26
	104(9)-201		100
Tissue SV MAXI (plus!)**	104(9)-310	spin/ vacuum	10
	104(9)-326		26
Blood SV mini	105-101	spin/ vacuum	100
	105-152		250
Blood SV Midi**	105-226	spin/ vacuum	26
	105-201		100
Blood SV MAXI**	105-310	spin/ vacuum	10
	105-326		26

Products	Cat. No.	Type	Size
GeneAll® Exgene™ for isolation of total DNA			
Cell SV mini	106-101	spin/ vacuum	100
	106-152		250
Cell SV MAXI**	106-310	spin/ vacuum	10
	106-326		26
Clinic SV mini	108-101	spin/ vacuum	100
	108-152		250
Clinic SV MAXI**	108-310	spin/ vacuum	10
	108-326		26
Plant SV mini	117-101	spin/ vacuum	100
	117-152		250
Plant SV Midi**	117-226	spin/ vacuum	26
	117-201		100
Plant SV MAXI**	117-310	spin/ vacuum	10
	117-326		26
GMO SV mini	107-150	spin/ vacuum	50
	107-102		200
GMO SV MAXI**	107-310	spin/ vacuum	10
	107-326		26
GeneAll® GenEx™ for isolation of total DNA			
Genomic Sx [†]	208-001	mini/ solution	100
Genomic Sx [†]	208-005	mini/ solution	500
Genomic Lx ^{††}	208-301	MAXI/ solution	100
GeneAll® RiboEx™ for preparation of total RNA			
RiboEx	301-001	solution	100
	301-002		200
RiboEx_column	301-150	Spin	50
Ribo_clear	303-150	Spin	50
GeneAll® Enzyme			
Taq DNA polymerase	501-025		250 U
	501-050		500 U
	501-100		1000 U
α-Taq DNA polymerase	502-025		250 U
	502-050		500 U
	502-100		1000 U
Pfu DNA polymerase	503-025		250 U
	503-050		500 U
	503-100		1000 U

* GeneAll® Tissue SV mini, Midi, and MAXI plus! kit provide the additional methods for the purification from animal whole blood.

** GeneAll® SV Midi/MAXI kits require the centrifuge which has a swinging-bucket rotor and ability of 4,000–5,000 xg.

[†] On the basis of DNA purification from 300 ul whole blood.

^{††} On the basis of DNA purification from 10 ml whole blood.

GeneAll® Expin™ Gel SV Brief protocol

** Before use this protocol, we strongly recommend you read first carefully the detailed protocol at page 12.*

Prepare a water bath or a heating block to 50°C

Check if precipitates appear in buffer GB. If so, dissolve completely at 50°C before proceeding.

1. Excise the DNA band of interest as small and rapid as possible.
2. Add 3 vol (ul) of buffer GB to 1 vol (mg) of gel.
3. Incubate at 50°C with occasional mixing until the gel is completely melted.
⇒ Check that the color of mixture is yellow after melting. If not, see the step 4 in detailed protocol at page 13.
4. Transfer the mixture to a SV column, centrifuge for 1 min, discard the filtrate.
5. Add 700 ul of buffer NW to the SV column, centrifuge for 30 sec and discard the filtrate.
6. Centrifuge for an additional 1 min and transfer the SV column to a new 1.5 ml tube.
7. Apply 50 ul of buffer EB or ddH₂O, let stand for 1 min and centrifuge for 1 min.

GeneAll® Expin™ PCR SV Brief protocol

** Before use this protocol, we strongly recommend you read first carefully the detailed protocol at page 12.*

1. Add 5 vol of buffer PB to 1 vol of the sample and mix thoroughly.
2. Transfer the mixture to a SV column, centrifuge for 30 sec and discard the filtrate.
3. Add 700 ul of buffer NW to the SV column, centrifuge for 30 sec and discard the filtrate.
4. Centrifuge for an additional 1 min and transfer the SV column to a new 1.5 ml tube.
5. Apply 50 ul of buffer EB or ddH₂O, let stand for 1 min and centrifuge for 1 min.

GeneAll® Expin™ CleanUp SV Brief protocol

** Before use this protocol, we strongly recommend you read first carefully the detailed protocol at page 12.*

Check if precipitates appear in buffer NR. If so, dissolve completely at 50°C before proceeding.

1. Add 10 (5) vol of buffer NR to 1 vol of the sample and mix thoroughly.
⇒ If the length of target DNA is larger than 100 bp, use 5 vol of buffer NR.
2. Transfer the mixture to a SV column, centrifuge for 30 sec and discard the filtrate.
3. Add 700 ul of buffer NW to the SV column, centrifuge for 30 sec and discard the filtrate.
4. Centrifuge for an additional 1 min and transfer the SV column to a new 1.5 ml tube.
5. Apply 50 ul of buffer EB or ddH₂O, let stand for 1 min and centrifuge for 1 min.



neAll

www.geneall.com

Banseok Bld., 128 Oguem-dong,
Songpa-gu, Seoul, KOREA 138-859

E-MAIL sales@geneall.com

T E L 82-2-407-0096

F A X 82-2-407-0779

©2009 GeneAll Biotechnology, All right reserved