

**Tribo™ Chromatin Immunoprecipitation (ChIP) Assay (Catalog# TBS8050)****DESCRIPTION**

**Chromatin Immunoprecipitation (ChIP)** is a powerful tool for analysis of interaction of proteins associated with chromosomal DNA. These proteins can be histone subunits and post-translational modification, or other chromatin associated proteins such as transcription factors, chromatin regulators. Additionally, ChIP can be used to identify regions of the genome associated with these proteins, or conversely, to identify proteins associated with a particular region of the genome.

**Tribo™ ChIP Assay Kit** is designed to use analysis of 25 assays for chromatin modification in tissue or cell.

**KIT CONTENTS****ChIP Kit:**

Protein A/G Magnetic Beads: 500 µL

Cell Lysis Buffer: 15 mL

Nuclear Lysis Buffer 15 mL

ChIP Dilution Buffer: 25 mL

Wash Buffer 1 (Low Salt): 25 mL

Wash Buffer 2 (High Salt): 25 mL

Wash Buffer 3 (LiCl): 25 mL

Wash Buffer 4 (TE): 25 mL

ChIP Elution Buffer: 1.5 mL

10x Glycine (1.25M), 5mL

Proteinase K (20 mg/mL): 0.5 mL (4°C)

Protease inhibitor Cocktail II (PI): 0.2 mL (-20 °C)

**Chromatin Magnetic Purification Kit**

Magnetic Binding Buffer A: 10.0 mL

Washing Buffer B: 7.0 mL (*Note: Add 28 mL of Ethanol before first-time use*)

Elution Buffer C: 2.0 mL

**OTHER COMPONENTS REQUIRED BUT NOT INCLUDED AS PART OF KIT**

Antibodies for ChIP, and negative control.

1 x PBS

Ethanol

**APPLICATIONS**

Chromatin modification for tissue or cell.

**STORAGE CONDITIONS**

The Kit stored at 4°C, and stable for 6 months from date of shipment.

**PROTOCOL****A: Tissue Process**

1. Quickly mince tissue using razor blades in 1.0 mL cold-PBS-PI in 10-cm culture dish on ice. Transfer the minced tissue to a new 2 ml microtube.
2. Mash the tissue as much as possible with the 1-ml pipette.
3. Simply spin the sample at 4°C and discard the supernatant and keep the pellet in the tube.

**B: X-linking the chromatin**

4. Add 1ml 1% formaldehyde in PBS-PI in each tube and resuspend the pellet. Incubate the sample for 10 min at room temperature. (1% Formaldehyde-PI: add 135µL of 37% formaldehyde to 5 mL of 1xPBS-PI).
5. Stop x-linking by adding 100 µL of 10 x glycine to the sample (10 times dilution of 1.25M glycine to final concentration of 0.125M). Incubate for 5 min at room temperature.
6. Simply spin the linked sample and pour off the supernatant and keep the pellet.
7. Wash 2x times with cold PBS-PI (1 ml). (The tissue may be snap frozen at this stage in store -80°C.)

**C: Lysis and Sonication**

8. Add 0.5 ml of Cell Lysis Buffer with PI to each pellet. Resuspend by mild pulses on the vortex mixer. Incubate on ice for 15 min; vortex every 5 min.
9. Simply spin at 2000x g, 4°C for 2min
10. Remove the supernatant and resuspend the cell in 0.5 mL nuclei lysis buffer (*Note: warm up nuclei lysis buffer to dissolve SDS precipitation*).
11. Homogenize the sample with glass Dounce Homogenizer (1ml) on ice for 10 strokes. Incubate the tube on ice for at least 10 min.
12. Sonicate the sample on ice to shear DNA to an average fragment size of 200~ 1000 bp. This will need optimizing as different equipment and tissue.
13. Centrifuge the sonicated chromatin at 10,000xg at 4 °C for 10 min (Do not exceed 15000 x g to prevent loss of chromatin). The supernatant is chromatin of the sample.
14. Remove supernatant to fresh microfuge tubes in 50- to 200 µL aliquots. Each 50 uL aliquot contains 1x 10<sup>6</sup> cell equivalents of lysate which is enough for one immunoprecipitation. Proceed to the next step or store at -80°C for 2 months.

15. To determine the DNA concentration, transfer 5 µl of the purified DNA into a tube containing 995 µl TE to give a 200-fold dilution and read the OD260. The concentration of DNA in µg/mL is OD260 x10,000. This used to calculate the DNA concentration of the chromatin preparation.
16. Perform agarose gel analysis of sonication (optional) Remove a 10-25 µL aliquote, and incubate for 10 min at 95 C. Add 1 µL Proteinase K and incubate at 62C for 2 hours with shaking. Load 10 µL or 20 µL on a 2% agarose gen with a 100-bp DNA marker.

**D. Immunoprecipitation**

17. Prepare enough ChIP Dilution Buffer containing protease inhibitors for the number of desired immunoprecipitations and store on ice. Add dilution buffer to each sample as shown below:

Sample Vol.	Dilution Buffer Vol.	PI Vol.
50 µL	450 µL	2.25 µL
100 µL	900 µL	4.5 µL

18. ChIP reaction includes the negative control and the antibody of interest (user supplied). It is recommended that the negative control IgG be from the same species as the antibody of interest.
19. Remove 5 µL of the supernatant as Input and save at 4°C until Section E (Reverse Cross-links of Protein/DAN Complexes to Free DNA).
20. Add target antibody (5 µg in 500 µL total volume) to each reaction tube.
21. For the negative control: add 1.0 µg of negative control antibody per tube.
22. Add 20 µL of fully resuspended protein A/G magnetic Beads to each reaction tube (Note: *Magnetic bead Slurry is well mixed before removing appropriate volume for IP, as magnetic bead will settle on the bottom of the tube over a short period of time*).
23. Incubate for 1 hour to overnight at 4°C with rotation.
24. Pellet Protein A/G magnetic beads with the magnetic separator or spin at 1500 rpm for 30 seconds. Remove the supernatant completely.
25. Wash the protein A/G bead-antibody/ chromatin complex by resuspending the beads in 0.5 mL each of the cold buffers in the order listed below and incubate for 3-5 min on a rotating platform between each wash followed by magnetic separation and careful removal of the supernatant fraction:  
*Low Salt Wash Buffer 1, one wash.*  
*High Salt Wash Buffer 2, one wash*  
*LiCl Wash Buffer 3, one wash.*  
*TE Buffer 4, one wash*

**E. Isolate DNA**

26. Warm the ChIP Elution Buffer to room temperature to ensure the SDS is dissolved in solution before proceeding. Prepare the final elution buffer for all samples including input (prepared in Section D, Step 18). For each tube, prepare elution buffer as follows:

ChIP Elution Buffer	100 µL
Proteinase K	1 µL

27. Incubate at 62 °C for 2 hours with shaking. Then, incubate at 95 °C for 10 min. Cool the samples down to room temperature.
28. Separate beads using a magnet separation device or spin at 2000 rpm for 30 seconds to drop down the beads to bottom of the tube. Carefully transfer the supernatant to a new tube.

**F. DNA Purification Using Magnetic Beads**

29. Add 0.3 mL of Magnetic Bind Reagent A to each 100 µL DNA sample tube, and mix, incubate at RT for 5-10min (**Note: 3 volume s of Bind Reagent A is used for every 1 volume of sample**).
30. Separate beads using a magnet separation device or spin at 2000 rpm for 30 seconds to drop down the beads to bottom of the tube. Carefully remove the supernatant (*Note: Do not disturb the beads*).
31. Wash the magnetic beads for 2 times with 500 µl of Wash Reagent B. Separate beads using a magnet separation device or spin at 2000 rpm for 30 seconds to drop down the beads to bottom of the tube. Carefully remove the supernatant (*Note: Do not disturb the beads*).
32. Dry the beads on the magnetic separation device at RT for 10 – 15 minutes.
33. Add 25-50 µL Elution Buffer C or DEPC water, and gently resuspend the magnetic particles.
34. Incubate the tube at 60°C for 5 minutes.
35. Separate beads using a magnet separation device or spin at 2000 rpm for 30 seconds to drop down the beads to bottom of the tube.
36. The eluate is purified DNA. It can be analyzed immediately or stored frozen at -20 °C. The DNA can be used for PCR analysis.

**RELATED PRODUCTS**

2xSyber Green Super qPCR Mix (Catalog# TBS4001)  
 2x Taqman Probe PCR Mix Catalog# TBS4002)

**This product is for *in vitro* research use only and is not intended for use in humans or animals in therapeutic and diagnosis purpose.**