

2x Fast Taqman Probe qPCR Master Mix

Fast Activation within 1 min

DESCRIPTION

Tribo™ 2x Fast Taqman Probe qPCR Master Mix with different ROX level and package size provides a convenient choice suitable for different qPCR instruments. It is a 2x concentrated premix containing Taq DNA polymerase, dNTPs, Mg²⁺, ROX reference dye and PCR buffer components. The proprietary chemical modification of the enzyme makes [PCR hot start activation much faster \(less than 1 min\)](#), do not need 10 min activation in the standard real-time qPCR amplification, which confers a significant reduction in non-specific PCR amplification ([please see the diagram in the data sheet](#)). As the enzyme buffer system is optimized for real time PCR, the product offers high amplification efficiency and high detection sensitivity in real-time PCR.

KEY FEATURES

- ❖ Fast: [Only 1min](#), Not need 10 or 5 min for Enzyme activation.
- ❖ Multiple choices: 2X Master Mix with different ROX level can be suitable for different instruments.

CATALOG AND KIT SIZE

STORE AT -20 °C

Kit Catalog	Description	Kit Size	qPCR Instrument
TBS4002R-1.25	2x Fast Taqman Probe qPCR Master Mix ROX	1x1.25 mL	ABI® PRISM® 7000, 7700, 7900HT, ABI® 7300 qPCR Systems, GeneAmp® 5700, StepOne™, and the StepOnePlus™.
TBS4002R-2.5	2x Fast Taqman Probe qPCR Master Mix ROX	2x1.25 mL	
TBS4002R-5	2x Fast Taqman Probe qPCR Master Mix ROX	4x1.25 mL	
TBS4002R-10	2x Fast Taqman Probe qPCR Master Mix ROX	1x10 mL	
TBS4002R-20	2x Fast Taqman Probe qPCR Master Mix ROX	2x10 mL	

Kit Catalog	Description	Size	qPCR Instrument
TBS4002LR-1.25	2x Fast Taqman Probe qPCR Master Mix Low ROX	1x1.25 mL	ABI® 7500 qPCR Systems, ViiA™ 7, QuantStudio™ 12K Flex, Agilent Mx3000P™ Mx3005P™ and Mx4000™
TBS4002LR-2.5	2x Fast Taqman Probe qPCR Master Mix Low ROX	2x1.25 mL	
TBS4002LR-5	2x Fast Taqman Probe qPCR Master Mix Low ROX	4x1.25 mL	
TBS4002LR-10	2x Fast Taqman Probe qPCR Master Mix Low ROX	1x10 mL	
TBS4002LR-20	2x Fast Taqman Probe qPCR Master Mix Low ROX	2x10 mL	

Kit Catalog	Description	Size	qPCR Instrument
TBS4002NR-1.25	2x Fast Taqman Probe qPCR Master Mix No ROX	1x1.25 mL	BioRad iCycler MiniOpticon, Opticon 2, Chromo 4, iQ5; Roche LightCycler 480; MJ Research DNA Engine Opticon 2, Chromo 4; Corbett Rotogene 3000, 6000
TBS4002NR-2.5	2x Fast Taqman Probe qPCR Master Mix No ROX	2x1.25 mL	
TBS4002NR-5	2x Fast Taqman Probe qPCR Master Mix No ROX	4x1.25 mL	
TBS4002NR-10	2x Fast Taqman Probe qPCR Master Mix No ROX	1x10 mL	
TBS4002NR-20	2x Fast Taqman Probe qPCR Master Mix No ROX	2x10 mL	

PROTOCOL

1. Set up the reaction Mixture in 20 µL reaction volume:

Components	Volume	Final Concentration
2xFast Taqman Probe qPCR master Mix	10 µL	1x
Forward Primer	Variable	100-600 nM
Reverse Primer	Variable	100-600 nM
Taqman Probe	Variable	50-300 nM
Template DNA	Variable	≤ 10 ng /reaction
Nuclease-Free H ₂ O	To 20 µL	

2. PCR amplification condition:

95°C, 1min → (95°C, 10S→60°C, 30S.) for 40 cycles.

STORAGE CONDITIONS

The kit is shipped on ice, and should be stored at -20°C for long-term storage. Shelf life of 12 months after receipt.

RELATED PRODUCTS

- 2x Fast Sybr Green qPCR Master Mix-ROX (TBS4001R)
- 2x Fast Sybr Green qPCR Master Mix-Low ROX (TBS4001LR)
- 2x Sybr Green Probe qPCR Master Mix-No ROX (TBS4001NR)
- All-in-One Reverse Transcription Reaction Kit (TBS4006)

REFERENCES


1. Birch DE. (1996): Simplified hot start PCR Nature, 381(6581): 445-446
2. Kellogg DE et al (1994): TaqStart Antibody: "hot start" PCR facilitated by a neutralizing monoclonal antibody directed against Taq DNA polymerase. Biotechniques. 16(6): 1134-7

Web site:

www.tribioscience.com

Research use only

2x Fast Taqman Probe qPCR Master Mix Amplification



2x Standard Taqman qPCR Master Mix Amplification

