

### INTRODUCTION

Tribioscience's Mitochondrial Membrane Potential Assay is based upon JC-10 dye, which is a superior alternative to JC-1 dye. JC-10 dye is a cationic, lipophilic dye that gets concentrated in healthy mitochondria to form reversible red-fluorescent JC-10 aggregates (Ex/Em=540/590nm). In apoptotic cells, the mitochondrial membrane potential collapses, which results in failure to retain JC-10 dye and the dye returns to its monomeric green form (EX/Em= 490/525 nm) in the cytosol.

The kit, not only measures mitochondrial membrane potential, but can be used for monitoring apoptosis and for screening of apoptotic activators and inhibitors.

### KIT CONTENT AND STORAGE CONDITIONS for 100 Assays

PART	PART#	Volume
JC-10 Dye Solution [200x]	TBS2049A	100 µl
MMP Assay Buffer [5 X]	TBS2049C	10 mL

Store the unopened kit at -20 °C stable for 12 months.  
Do not use past kit expiration date.

### PRECAUTIONS

Wear protective gloves, clothing, eye, and face protection. Wash hands thoroughly after handling.

### REAGENT PREPARATION

**Bring all reagents to room temperature before use.**

- **MMP Assay Buffer:** Dilute the MMP-Assay Buffer [5X] to 1X with molecular grade water in ratio 1:4 to get a desired volume of MMP-Assay Buffer [1X] for the assay.
- **JC-10 Dye Working Solution:** Centrifuge the JC-10 Dye vial briefly. Prepare the desired volume of JC-10 Dye Working Solution by adding dye to prewarmed cell culture medium in ratio 1:100. For example, add 10 µl of JC-10 dye (100x) to 1 ml of prewarmed cell culture medium.

### ASSAY PROCEDURE

1. Culture cells in 96-well black microtiter plate at a density of  $5 \times 10^4$  -  $5 \times 10^5$  cells per well in 100 µl cell culture medium overnight in incubator (5% CO<sub>2</sub>, 37°C).
2. Treat the cells with or without the test compound in duplicate manner.
3. Incubate the plates in incubator (5% CO<sub>2</sub>, 37°C) for a desired period time, like 4-6 hours. NOTE: The optimum incubation time may vary depending upon cell type and compound property.
4. Replace the cell culture medium with 100 µl of JC-10 Dye Working Solution under cell culture hood.
5. Incubate the plates in incubator (5% CO<sub>2</sub>, 37°C) for 15-60 minutes. NOTE: The optimum incubation time may vary depending upon cell type and cell concentration.
6. Aspirate the supernatant and add 100 µl of prewarmed 1 X MMP-Assay Buffer per well.
7. Read the fluorescence at Ex/Em: 535nm/595 nm for JC-10 aggregates (red fluorescence) and at Ex/Em: 485 nm/535 nm for JC-10 monomers (green Fluorescence).
8. Plot the graph against the test compound concentration and ratio of red to green fluorescence to determine the effects of compound on cell health.

**Research Use only.**