Aspergillus Genomic DNA Isolation Kit Catalog# TBS6009

DESCRIPTION

The Aspergillus Genomic DNA Isolation Kit is designed for reliable and efficient isolation of genomic DNA from a broad spectrum of fungus species, including Aspergills fumigatus, Aspergills nidulans, Aspergills nivens var. aureus, Candidaalbicans, Pichia pastoris, Saccharomyces cerevisiae, Schizo- saccharomyces pombe, and any fungi whose cell walls are susceptible to yeast lytic enzyme lysis. The resultinggenomic DNA can be used directly for all molecular biology analysis such as Southern Blot, PCR, and restriction enzyme digestion.

KIT COMPONENTS (100 TESTS)

Part#	Name	Size	Store
B1	Extraction Buffer 1	50 mL	RT
B2	Extraction Buffer 2	1 mL	-20 °C
B3	Digest Buffer 3	0.5 mL	-20 °C
B4	Purification Buffer 4	25 mL	4 °C
B5	TE buffer	2 mL	4 °C
Shelf-life: 6 month after receipt.			

APPLICATIONS

- PCR.
- Southern Blot.
- Restriction Enzyme Digestion

REAGENTS & EQUIPMENT NOT PROVIDED

- Isopropanol
- Ethanol

SAMPLE TYPES

• Aspergillus

DNA ISOLATION PROCEDURE

- **1.** 0.1- 0.2 g of fresh fungal mycelium was transferred to a sterilized 2 ml Eppendorf tube with 0.5 ml of Extraction Buffer 1.
- **2.** Add 10 μL Extraction Buffer 2. Mix well. Incubatefor 30 min. at 37°C.
- **3.** Add 5 μL Digest Buffer 3. Mix well. Placed in awater bath at 65°C for 0.5-1 hr.
- **4.** Then place the sample to centrifuged at 10,000 rpmfor 10 min at room temperature;
- **5.** Collect the supernatant into a new microtube (2mL).
- 6. Add equal volume of Purification Buffer 4, then mixvery well. Centrifuge again at 10,000 rpm for 5 min at RT;
- 7. Collect the supernatant into a new tube (2 mL), andthen, add 1 volume of ice-cold isopropanol, and incubated at -20°C for 1h or -80°C for 30 min; Centrifuge for 15 min at 4°C, 13,000 rpm;
 (PS: The equal volume of ice-cold isopropanol canreplace by 2 volume of

isopropanol canreplace by 2 volume of ethanol.)

- **8.** Supernatant is decanted, and DNA pellet is washedwith 1ml of 70% ethanol, repeat 3 times;
- **9.** DNA pellet was air-dried and dissolved in 100 µLdistilled water or TE buffer.

RELATED PRODUCTS

Universal Aspergillus qPCR Detection Kit (TBS42020) Aspergillus Flavus qPCR Detection Kit (TBS42021) Aspergillus Fumigatus qPCR Detection Kit (TBS42022) Aspergillus Niger qPCR Detection Kit (TBS42023) Aspergillus Terreus qPCR Detection Kit (TBS42024)