


7-MINUTE DNA ISOLATION FROM MOUSE TAIL AND OTHER TISSUES

TBS6005

Mouse tail DNA extraction kit provides a simple and rapid method to isolate genomic DNA from mouse tail for genotyping. The formulated buffer system is a ready-to-use single reagent without other preparation. No organic extraction and alcohol precipitation are needed and multiple samples can be easily processed simultaneously. It incubates with sample just for **7 minutes** instead of overnight incubation.

Features of DNA Isolation Kit:

- Pre-mixed buffer system
- Ready for DNA in **just 7 minutes**
- Simple Procedure: Add-incubate
- Stable and consistent result
- No need of additional reagents
- Optimized protocols



Company: TriBioScience Inc; Address: 365 San Aleso Ave, Sunnyvale, CA 94085

Phone: 408-498-0197/650-917-9269; Fax: 650-618-5498

Email: info@tribioscience.com; Web: www: tribioscience.com