



# Fast Mouse Genotyping System (TBS4033)


Tel: 408-498-0197 | E-mail: [order@tribioscience.com](mailto:order@tribioscience.com) | Fax: 650-618-5498  
 365 San Aleso Ave. Sunnyvale, CA 94085; [www.tribioscience.com](http://www.tribioscience.com)

Fast Mouse Genotyping System is designed for rapid extraction and PCR amplification of genomic DNA from mouse tail/toe/ear punch and other tissues. It contains all reagents required for DNA extraction and PCR amplification. PCR products can be directly loaded onto agarose gels to visualize amplified target.

## Features

1. Fast DNA extraction only 10 min. No overnight
2. Multiplexing: strip tube for parallel processing
3. Cost-Saving:  
 Less than \$0.79/sample;  
 Sigma: \$2.63/sample.
4. PCR products can be directly analyzed in gel
5. Better than brand products like sigma.

### One Step DNA Isolation in 10 minutes



### Simple PCR Amplification ( 20 µL/ Reaction)





Figure 1. Sigma PCR Mix



Figure 2. Tribio\_2xGenotyping PCR

### Product Package for 250 Samples

Mouse Tail Fast DNA Extraction

2x Genotyping PCR Ready Mix



For more information, please go to [www.tribioscience.com](http://www.tribioscience.com)

Follow us on Facebook

<https://www.facebook.com/tribioscience/>

**Fast: Fidelity, Accurate, Simple and Time-saving**