

## 10x Tris Buffered Saline (10x TBS, TBS5024)

---

### DESCRIPTION

Tris buffered Saline (TBS) is commonly used as a wash buffer and diluent solution for Western blot, Elisa, and general wash buffer in immunoassay applications.

Some applications need the detergent Tween-20 in the buffer. So, we provide several different formulations with or without Tween-20 in the buffer to make an easy choice for the application.

The buffer concentration and ingredients are listed below:

#### **10X Formulation (TBS5024):**

1500 mM Sodium Chloride

500 mM Tris-HCl

Supplied at pH 7.2 - 7.4

### DIRECTIONS FOR USE

Simple dilute the concentrated TBS (10x) to a 1X solution using ddH<sub>2</sub>O.

### STORAGE

Store at room temperature.

**Orders** Tel:408-498-0197  
order@tribioscience.com

**Support** Tel:408-498-0197  
info@tribioscience.com

**Web:** www.tribioscience.com

**For research use only.**