
DEPC-Treated H₂O

DESCRIPTION

DEPC-treated water is suitable for use with RNA. It is prepared by incubating RNase and DNase-free water with 0.1% diethylpyrocarbonate (DEPC). The DEPC is then removed in the autoclave.

Sterile filtered. No DNase or RNase activity detected. Molecular Biology Grade.

STORAGE

Supplied as ready-to-use RNase and DNase-Free H₂O, 500mL. Shipping and storage at room temperature.

APPLICATIONS

This item is used for molecular biology assay in the RNase and DNase-free conditions.

RELATED PRODUCTS

6x Protein Loading Buffer (Catalog# TBS5014)

10x PBS (Catalog# TBS5027)

Stripping Buffer (Catalog# TBS5020)

RIPA Buffer (Catalog# TBS5017)

Protein Assay Kit (Catalog# TBS2005)

2x qPCR Hot Start Super Mix (Catalog# TBS4001)

Tribo™ Taq DNA Polymerase (Catalog# TBS4010)

This product is for *in vitro* research use only and is not intended for use in humans or animals in therapeutic or diagnostic procedures.

Contact

Phone: 408-498-0197

Email: info@tribioscience.com

Fax: 650-618-5498

Web: www.tribioscience.com