

DESCRIPTION

Tris buffered Saline (TBS) is commonly used as a wash buffer and diluent solution for Western blot, Elisa, and general wash buffer in immunoassay applications.

Some applications need the detergent Tween-20 in the buffer. So, we provide several different formulations with or without Tween-20 in the buffer to make an easy choice for application.

The buffer concentration and ingredients are listed below:

10x TBS Formulation with Tween-20 (TBS5008T):

1500 mM Sodium Chloride
500 mM Tris-HCl
0.5% Tween-20
Supplied at pH 7.2 - 7.4

DIRECTIONS FOR USE

Simple dilute concentrated TBST(10x) to a 1X solution using ddH₂O.

STORAGE

Store at room temperature.

Orders

Tel:408-498-0197
order@tribioscience.com

Support

Tel:408-498-0197
info@tribioscience.com

Web

www.tribioscience.com

For research use only.