

20% Sucrose Gradient Solution in STE (pH 7.4)

Catalog
TBS5075

Unit Size
500 mL

DESCRIPTION

Tribioscience's 20% Sucrose Gradient Solution is ready to use for nucleic acid (RNA and DNA), viral vector purification in gradient centrifuge.

We can provide customized specific products to meet customer needs.

COMPOSITION

20% Sucrose

50 mM Tris-HCl

100 mM NaCl

0.5 mM EDTA

pH Value: 7.4 (7.3 ~ 7.5)

Sterile filtered with 0.2um filter.

STORAGE

Store at 2-8°C. Use for 12 months from date of receipt.

APPLICATION

- Calcium transfection
- Calcium Phosphate Transfection
- Mammalian cell transfection
- Virus production
- Lentivirus and retrovirus transfection
- AAV transfection

Research use only.