

**Tris Buffered Magnesium Acetate Saline (TMS, 1x); Catalog: TBS5085**

---

**DESCRIPTION**

Tris buffered Magnesium Acetate Saline (TMS) is widely used as a wash buffer or diluent solution for Western blot, Elisa, and general wash buffer in immunoassay applications.

The Buffer is a ready- to-use format without any dilution.

The buffer concentration and ingredients are listed below:

**1X Formulation (TBS5085):**

150 mM Sodium Chloride.

50 mM Tris-acetate.

50 mM Magnesium acetate.

Supplied at pH 7.2 - 7.6

**Size:** 1L.

**STORAGE**

Store at room temperature.

**Orders** Tel:408-498-0197

order@tribioscience.com

**Support** Tel:408-498-0197

info@tribioscience.com

**Web:** www.tribioscience.com

For research use only.