

DESCRIPTION

Tris buffered Magnesium Acetate Saline with Tween-20 (TMST) is widely used as a wash buffer or diluent solution for Western blot, Elisa, and general wash buffer in immunoassay applications.

The Buffer is a ready- to-use format without any dilution.

The buffer concentration and ingredients are listed below:

1X Formulation (TBS5085):

150 mM Sodium Chloride.

50 mM Tris-acetate.

50 mM Magnesium acetate.

0.05% Tween-20

Supplied at pH 7.2 - 7.6

Size: 1L.

STORAGE

Store at room temperature.

Orders Tel:408-498-0197

order@tribioscience.com

Support Tel:408-498-0197

info@tribioscience.com

Web: www.tribioscience.com

For research use only.