
DEPC-Treated H2O

DESCRIPTION

DEPC-treated water is suitable for use with RNA. It is prepared by incubation with 0.1% diethyl pyrocarbonate (DEPC) and then autoclaved to remove the DEPC. Sterile filtered. No DNase, RNase, or protease activity detected. Molecular Biology Grade.

STORAGE

Supplied as ready-to-use RNase and DNase-Free H₂O, 500 mL. Shipping and store at RT.

APPLICATIONS

This item is used for molecular biology assays under RNase-, DNase-, and protease-free conditions.

RELATED PRODUCTS

6x Protein Loading Buffer (catalog# TBS5014)
10x PBS (catalog# TBS5027)
Stripping Buffer (catalog# TBS5020)
RIPA Buffer (catalog# TBS5017)
Protein Assay Kit (catalog# TBS2005)
2x qPCR Hot Start Super Mix (catalog# TBS4001)
Tribo™ Taq DNA Polymerase (catalog# TBS4010)

This product is for *in vitro* research use only and is not intended for use in humans or animals in therapeutic or diagnostic procedures.

Contact

Phone: 408-498-0197

Email: info@tribioscience.com

Web: www.tribioscience.com