

Catalog
TBS5066-1L

Unit Size
1 L

Description

Tribio Cleaning Solution is designed for the BD Accuri C6 Flow Cytometer. It effectively removes residual cells, proteins, and fluorescent dyes from the flow cell and fluidic system, ensuring optimal performance and accurate results.

Application

- Routine maintenance and cleaning of BD Accuri C6 Flow Cytometer.
- Removal of sample residues after high-concentration or sticky samples.
- Prevention of clogs and contamination in the fluidic system.

Main Features

- Removes residual cells, proteins, and fluorescent dyes.
- Maintains stable fluidic flow and accurate readings.
- Extends instrument life by preventing blockages and pump strain.
- Compatible with BD Accuri C6 system.
- Safe for routine use and non-corrosive to instrument components.

Storage

Store at room temperature (20-25°C).

Relative Products

TBS2005	Protein Assay Kit
TBS5001	Protein Lysis Buffer
TBS5003-0	1xPBS
TBS5011-10	10 x PBST
TBS5008	TBS
TBS5008T	TBST
TBS5037	Urea Solution (8M)
TBS5038	Ammonium Sulfate Solution
TBS5040	0.5M EDTA Solution
TBS5041	MOPS Buffer
TBS5048	2% BSA in PBS
TBS5049-2	2% BSA in TBST
TBS5049-5	5% BSA in TBST
TBS5057	0.1% BSA ELISA Assay Buffer
TBS5077	Antigen retrieval Citrate Solution-pH6.0
TBS8024	MCF10 Cell Complete Medium
TBS8027	EMEM with L-Glutamine
TBS8030	McCoy's 5A Medium
TBS8031	30% BSA Solution in DPBS
TBS8047	MCF-7 Cell Complete Medium
TBS8056	LB Medium
TBS8057	SOB Medium
TBS8058	SOC Medium
TBS8059	2x YT Medium
TBS8060	Terrific Broth (TB) Medium
TBS8061	DMEM-high Glucose Medium

This product is for research use only.