

Catalog
TBS5067-1L

Unit Size
1 L

Description

Tribio Decontamination Solution is specially formulated for the Accuri C6 Flow Cytometer to remove residual cells, proteins, and microbial contaminants from the instrument's fluidics system. It ensures reliable operation, prevents cross-contamination, and maintains consistent data quality.

Application

- Cleaning and decontaminating the internal fluidics of Accuri C6 before or after experiments.
- Preventing cross-contamination between different samples.
- Maintaining optimal instrument performance and prolonging its operational life.

Main Features

- Effectively removes cellular debris, proteins, and microbial contaminants.
- Compatible with Accuri C6 sheath fluid and standard cleaning procedures.
- Help preserve accurate flow cytometry measurements.
- Non-corrosive to instrument components.
- Ready-to-use formulation, requiring no dilution.

Storage

Store at room temperature (20-25°C).

Relative Products

TBS2005	Protein Assay Kit
TBS5001	Protein Lysis Buffer
TBS5003-0	1xPBS
TBS5011-10	10 x PBST
TBS5008	TBS
TBS5008T	TBST
TBS5037	Urea Solution (8M)
TBS5038	Ammonium Sulfate Solution
TBS5040	0.5M EDTA Solution
TBS5041	MOPS Buffer
TBS5048	2% BSA in PBS
TBS5049-2	2% BSA in TBST
TBS5049-5	5% BSA in TBST
TBS5057	0.1% BSA ELISA Assay Buffer
TBS5077	Antigen retrieval Citrate Solution-pH6.0
TBS8024	MCF10 Cell Complete Medium
TBS8027	EMEM with L-Glutamine
TBS8030	McCoy's 5A Medium
TBS8031	30% BSA Solution in DPBS
TBS8047	MCF-7 Cell Complete Medium
TBS8056	LB Medium
TBS8057	SOB Medium
TBS8058	SOC Medium
TBS8059	2x YT Medium
TBS8060	Terrific Broth (TB) Medium
TBS8061	DMEM-high Glucose Medium

This product is for research use only.