




Fast Mouse Genotyping System (TBS4033)

Fast Mouse Genotyping System is designed for rapid extraction and PCR amplification of genomic DNA from mouse tail/toe/ear punch and other tissues. It contains all reagents required for DNA extraction and PCR amplification. PCR products can be directly loaded onto agarose gels to visualize amplified target.

Features

1. Fast DNA extraction only 12 min. No overnight
2. Multiplexing: strip tube for parallel processing
3. Cost-Saving:
Less than \$0.79/sample;
Sigma: \$2.63/sample.
4. PCR products can be directly analyzed in gel
5. Better than brand products like sigma.

One Step DNA Isolation in 12 minutes



Simple PCR Amplification (20 µL/ Reaction)





Figure 1. Sigma PCR Mix



Figure 2. Tribio_2xGenotyping PCR

Product Package for 250 Samples

- Mouse Tail Fast DNA Extraction
- 2x Genotyping PCR Ready Mix



For more information, please go to www.tribioscience.com

Follow us on Facebook

<https://www.facebook.com/tribioscience/>

Fast: Fidelity, Accurate, Simple and Itime-saving