

Tribo™ Human Insulin-like Growth Factor (Human IGF1) Primer-Probe Mix (Catalog# TBS4101)

DESCRIPTION

Tribo™ Human IGF1 Primer-Probe Mix is optimized to real-time PCR based on primer and probe. It contains primers (forward and reverse) and probe labeled with a FAM reporter dye at the 5' end of the probe. The mix can be used for detecting gene expression.

APPLICATIONS

This kit can be used for FAM dye-based real-time PCR.

KEY FEATURES

- Easy-to-use pre-mixed package
- High fidelity yield and sensitivity
- Robust and high throughput
- Work with difficult sequences

KIT CONTENTS

10x Primer-Probe Mix: 100 µL

It contains 18 uM of each primer and 5uM for the probe.

STORAGE CONDITIONS

The kit is shipped on ice, and should be stored at -20°C for long-term storage. Shelf life of 12 months after receipt.

AMPLICON SIZE

68 bp

RELATED PRODUCTS:

2x Taqman qPCR Super Mix (catalog# TBS4001)
Blood DNA Extraction kit (catalog# TBS6004)