

## Tribo™ Human Myostatin (Human MSTN) Primer-Probe Mix (Catalog# TBS4102)

### DESCRIPTION

Tribo™ Human MSTN Primer-Probe Mix is optimized to real-time PCR based on primer and probe. It contains primers (forward and reverse) and probe labeled with a FAM reporter dye at the 5' end of the probe. The mix can be used for detecting gene expression.

### APPLICATIONS

This kit can be used for FAM dye-based real-time PCR.

### KEY FEATURES

- Easy-to-use pre-mixed package
- High fidelity yield and sensitivity
- Robust and high throughput
- Work with difficult sequences

### KIT CONTENTS

10x Primer-Probe Mix: 100 µL

### STORAGE CONDITIONS

The kit is shipped on ice, and should be stored at -20°C for long-term storage. Shelf life of 12 months after receipt.

### AMPLICON SIZE

69 bp

### RELATED PRODUCTS:

2x Taqman qPCR Super Mix (catalog# TBS4001)  
Blood DNA Extraction kit (catalog# TBS6004)