

0.2M Zinc Chloride Solution (Catalog# TBS5047)**DESCRIPTION**

Zinc Chloride Solution is a formulated water solution. This ready-to-use solution can be used for protein structure and fold assays, nucleotide acid precipitation, and other biological analysis.

APPLICATIONS

Applications include protein structure and fold analysis, nucleotide acid precipitation, and other biological analysis.

CONCENTRATION AND SIZE

The solution contains **0.2 M Zinc Chloride**. It is supplied as **100 mL** in each bottle.

STORAGE CONDITIONS

Store at 4°C and protected from light. Shelf life: 12 months.

RELATED PRODUCTS

8M Urea Solution(#TBS5037)
4M Ammonium Sulfate Solution (#TBS5038)
1M DTT Solution (#TBS5039)
0.5M EDTA Solution (#TBS5040)
1M Sodium Citrate Solution (#TBS5043)
1M sodium Acetrate solution (#TBS5044)

This product is for *in vitro* research use only and is not intended for use in humans or animals in therapeutic or diagnostic procedures.

Contact

Phone: 650-917-9269

Email: info@tribioscience.com

Fax: 650-618-5498

Web: www.tribioscience.com