

## Tribio™ Water for Cell Culture (Catalog# TBS5050)

**DESCRIPTION**

Tribio™ Water for Cell Culture is high quality water for use as a solvent in the preparation of cell culture media and laboratory reagents.

Tribio™ Water for Cell Culture is prepared by distillation and filtration (0.1µm), with no added substances.

**APPLICATIONS**

Tribio™ Water for Cell Culture can be used for a wide variety of applications, including preparing powder or cell culture media, making solutions and diluting buffers.

**CONTENTS**

500 mL/each bottle.

**Sterility:** The water is sterilized and filtered with 0.1 µm filter.

**Storage conditions:** Store at room temperature. Shelf life: 24 months after receipt.

**This product is for *in vitro* research use only and is not intended for use in humans or animals.**

**This product is not intended for use as a therapeutic or in diagnostic procedures.**

**Please visit [www.tribioscience.com](http://www.tribioscience.com) for a complete listing of recommended companion products.**

**Orders** Tel: 650-917-9269  
[order@tribioscience.com](mailto:order@tribioscience.com)

**Support** Tel: 650-917-9269  
[info@tribioscience.com](mailto:info@tribioscience.com)

**RELATIVE PRODUCTS**

Germine Transmission Qualified Fetal Bovine Serum(FBS) (Catalog#: TBS8001)

ESC/iPSC-qualified Fetal Bovine Serum(FBS)( Catalog#: TBS8002)

0.1% Gelatin solution ( Catalog#: TBS8004)

PBS, 1x ( Catalog#: TBS5003)