

DMEM-high Glucose Medium (Catalog# TBS8061-500 mL)

With High glucose, L-glutamine and phenol red, without sodium pyruvate

DESCRIPTION

DMEM (Dulbecco's Modified Eagle Medium) is one of widely used basal mediums for supporting the growth of many different mammalian cells. DMEM is unique from other media as it contains 4 times the concentration of amino acids and vitamins than the original Eagle's Minimal Essential Medium. Cells successfully cultured in DMEM include primary fibroblasts, neurons, glial cells, HUVECs, and smooth muscle cells, as well as cell lines such as HeLa, 293, Cos-7, and PC-12. We offer a variety of DMEM modifications for a range of cell culture applications.

TBS8061-500ML is modified with 4,500 mg/L glucose, L-glutamine, phenol red, sodium bicarbonate without sodium pyruvate. Users can make any suitable modification based on the specific needs for different cell lines.

PACK SIZE

1x500 mL/bottle, Store at 2-8°C °C in dark till use.

RELATED PRODUCTS

MSC Medium (TBS8021)
 Chondrogenic Differentiation Medium (TBS8062)
 RPMI16040 Medium (TBS8063)
 ESC/iPSC-qualified FBS (TBS8002)
 Adipocyte Differentiation Cocktail (TBS8017)
 0.1% Gelatin Solution (TBS8004)
 1.25M Calcium Chloride (TBS5071)
 2.5M Calcium Chloride (TBS5072)
 2x HBS, pH7.05 (TBS5076)
 Cell Culture Grad Water (TBS5050)
 LB Medium (TBS8056)
 SOB Medium (TBS8057)
 SOC Broth Medium (TBS8058)
 2xYT Broth Medium (TBS8059)

Research use only.

COMPOSITIONS

Ingredients	Concentration (mg/L)
Glycine	30
L-Arginine hydrochloride	84
L-Cystine 2HCl	63
L-Glutamine	584
L-Histidine hydrochloride-H2O	42
L-Isoleucine	105
L-Leucine	105
L-Lysine hydrochloride	146
L-Methionine	30
L-Phenylalanine	66
L-Serine	42
L-Threonine	95
L-Tryptophan	16
L-Tyrosine disodium salt dihydrate	104
L-Valine	94
Choline chloride	4
D-Calcium pantothenate	4
Folic Acid	4
Niacinamide	4
Pyridoxine hydrochloride	4
Riboflavin	0.4
Thiamine hydrochloride	4
i-Inositol	7.2
Sodium Bicarbonate	3700
CaCl2 (anhyd.)	200
Fe(NO3)3·9H2O	0.1
MgSO4(anhyd.)	97.67
KCl	400
NaCl	6400
NaH2PO4-H2O	125
D-Glucose (Dextrose)	4500
Phenol Red	15