

## Calcium Chloride (CaCl<sub>2</sub>) Solution (1.25 M or 2.5M) For Transfection

---

Catalog	Concentration (M)	Unit Size
TBS5071	1.25	100 mL
TBS5072	2.5	100 mL

### DESCRIPTION

Tribioscience's CaCl<sub>2</sub> Solution at a concentration of 1.25M or 2.5M is used for chemical transfection. It can be used with 2xHBS, pH 7.05 TBS5076-7.05.

### COMPOSITION

CaCl<sub>2</sub> (1.25 M or 2.5M) in water.

Sterile filtered with 0.2 um filter.

### STORAGE

Room temperature (18-25°C). Use by 12 months from the date of receipt.

### APPLICATION

- Calcium transfection
- Calcium Phosphate Transfection
- Mammalian cell transfection
- Virus production
- Lentivirus and retrovirus transfection
- AAV transfection

**For research use only.**