

AntiMullerian Hormone (AMH) Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10013-0.25MG	0.25mg
TBS10013-1MG	1mg

Description: Mouse Monoclonal Antibody to Antimullerian Hormone (AMH).

Specificity: This antibody is specific to human AMH.

Host Animal: Mouse.

Isotype: IgG1.

Source: Mouse ascites.

Immunogen: Recombinant Human AMH.

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 4.52 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN₃, pH 7.2.

Pairing Assays: Pairs with and TBS10012.

Purity: >95% by HPLC and SDS-PAGE.

Application: For human AMH quantitative assays by EIA or CLIA and Lateral Flow.

Storage: 2-8°C for 1 year. Long-term store at -20°C.

For Research Use only