

## AntiMullerian Hormone (AMH) Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TMA055-0.25MG	0.25mg
TMA055-1MG	1mg

**Description:** Mouse Monoclonal Antibody to Antimullerian Hormone (AMH).

**Specificity:** This antibody is specific to human AMH.

**Host Animal:** Mouse.

**Isotype:** IgG1.

**Source:** Mouse ascites.

**Immunogen:** Recombinant Human AMH.

**Format:** Purified, liquid.

**Purification:** DEAE chromatography.

**Concentration:** 4.52 mg/ml.

**Affinity Constant:** Not determined.

**Buffer:** 10 mM PBS, 0.05% NaN<sub>3</sub>, pH 7.2.

**Pairing Assays:** Pairs with and TMA056.

**Purity:** >95% by HPLC and SDS-PAGE.

**Application:** For human AMH quantitative assays by EIA or CLIA and Lateral Flow.

**Storage:** 2-8°C for 1 year. Long-term store at -20°C.

**For Research Use only**