

Monoclonal Anti-Pregnancy Associated Plasma Protein A (PAPP-A) Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10173-0.25MG	0.25mg
TBS10173-1MG	1mg

Description: Mouse Monoclonal Anti-Pregnancy Associated Plasma Protein A (PAPP-A) Antibody.

Specificity: This antibody reacts specifically to human PAPP-A.

Host Animal:	Mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Human PAPP-A.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	4.97 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10174 for human PAPP-A detection.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human PAPP-A quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only