

## Monoclonal Anti-Pregnancy Associated Plasma Protein A (PAPP-A) Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TBS10174-0.25MG	0.25mg
TBS10174-1MG	1mg

**Description:** Mouse Monoclonal Anti-Pregnancy Associated Plasma Protein A (PAPP-A) Antibody.

**Specificity:** This antibody reacts specifically to human PAPP-A.

<b>Host Animal:</b>	Mouse.
<b>Isotype:</b>	IgG1.
<b>Source:</b>	Mouse ascites.
<b>Immunogen:</b>	Human PAPP-A.
<b>Format:</b>	Purified, liquid.
<b>Purification:</b>	DEAE chromatography.
<b>Concentration:</b>	4.35 mg/ml.
<b>Affinity Constant:</b>	Not determined.
<b>Buffer:</b>	10 mM PBS, 0.05% NaN <sub>3</sub> , pH 7.2.
<b>Pairing Assays:</b>	Pair with TBS10173 for human PAPP-A detection.
<b>Purity:</b>	>95% by HPLC and SDS-PAGE.
<b>Application:</b>	For human PAPP-A quantitative assays by EIA or CLIA and Lateral Flow.
<b>Storage:</b>	2-8°C for 1 year. Long-term store at -20°C.

**For Research Use only**