

Monoclonal Anti-Thyroid Stimulating Hormone (TSH) Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10183-0.25MG	0.25mg
TBS10183-1MG	1mg

Description:	Mouse monoclonal anti-human thyroid stimulating hormone (TSH) antibody.
Specificity:	TSH 100% hCG 0.0% Alpha-hCG 0.0% Beta-hCG 0.0% LH 0.1% FSH 0.1%.
Host Animal:	Mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Human TSH.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	4.95 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Paired with TBS10182.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human TSH quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only