

Monoclonal Anti-human Progesterone Antibody

Specification Sheet

Catalog Number	Unit Size
TMP104-0.25MG	0.25mg
TMP104-1MG	1mg

Description:	Mouse monoclonal anti-human progesterone antibody.
Specificity:	This antibody reacts specifically to human progesterone.
Host Animal:	Mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Progesterone-BSA.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	4.42 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	N/A.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human progesterone quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only