

Monoclonal Anti-human Cystatin C Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10090-0.2MG	0.2mg
TBS10090-1MG	1mg

Description: Monoclonal Anti-Cystatin C antibody. This item (TBS10090) can be paired with TBS10091 for ELISA.

Specificity:	This antibody reacts specifically with human cystatin C.
Host Animal:	Balb/c mouse.
Isotype:	IgG1, kappa light chain.
Source:	Mouse ascites.
Immunogen:	Human Cystatin C.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	4.36 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10091 for human cystatin C detection.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human cystatin C quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only