

## Monoclonal Anti-human Thyroglobulin (TG) Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TBS10185-0.2MG	0.2mg
TBS10185-1MG	1mg

**Description:** Monoclonal Anti-human thyroglobulin (TG) antibody. This item (TBS10185) can be paired with TBS10186 for ELISA.

<b>Specificity:</b>	This antibody reacts with human TG.
<b>Host Animal:</b>	Balb/c mouse.
<b>Isotype:</b>	IgG2a.
<b>Source:</b>	Mouse ascites.
<b>Immunogen:</b>	Human thyroglobulin.
<b>Format:</b>	Purified, liquid.
<b>Conjugate:</b>	No.
<b>Purification:</b>	DEAE chromatography.
<b>Concentration:</b>	6.9 mg/ml.
<b>Affinity Constant:</b>	Not determined.
<b>Buffer:</b>	10 mM Tris-HCL buffer, 0.85% NaCl, 0.1% NaN <sub>3</sub> , pH 8.0.
<b>Pairing Assays:</b>	Pair with TBS10186.
<b>Purity:</b>	>95% by HPLC and SDS-PAGE.
<b>Application:</b>	For human thyroglobulin quantitative assays by EIA or CLIA, Lateral Flow, and microfluidics.
<b>Storage:</b>	2-8°C for 1 year. Long-term store at -20°C.

**For Research Use only**