

## Monoclonal Anti-human Transferrin Antibody-HRP

### Specification Sheet

Catalog Number	Unit Size
TBS10190HRP-0.2MG	0.2mg
TBS10190HRP-1MG	1mg

**Description:** Monoclonal Anti-human transferrin antibody conjugated with HRP. This item (TBS10190HRP) as detective antibody can be paired with TBS10189 as capture antibody for ELISA. The ELISA performance is displayed in Figure 1.

**Specificity:** This antibody reacts specifically with human transferrin.

**Host Animal:** Balb/c mouse.

**Isotype:** IgG1.

**Source:** Mouse ascites.

**Immunogen:** Human transferrin

**Format:** Purified, liquid.

**Conjugate:** HRP.

**Purification:** Chromatography.

**Concentration:** 1.0 mg/ml.


**Affinity Constant:** Not determined.  $R^2=0.998$

**Buffer:** 10 mM PBS, 1% BSA, 50% glycerol, pH 7.2.

**Pairing Assays:** Pair with TBS10189 (capture antibody).

**Application:** For human transferrin quantitative assays by EIA, and microfluidics.

**Storage:** 2-8°C for 1 year.



**For Research Use only**