

Monoclonal Anti-human Ferritin Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10113-0.2MG	0.2mg
TBS10113-1MG	1mg

Description: Monoclonal Anti-human ferritin antibody. This item (TBS10113) as detective Mab can be paired with TBS10114 as capture Mab for ELISA.

Specificity:	This antibody reacts with human ferritin.
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Human liver ferritin.
Format:	Purified, liquid.
Conjugate:	No.
Purification:	DEAE chromatography.
Concentration:	5.72 mg/ml.
Affinity Constant:	Not determined.
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10114.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human ferritin quantitative assays by EIA or CLIA, Lateral Flow, and microfluidics.
Storage:	2-8°C for 1 year. Long-term store at -20°C.

For Research Use only