

## Monoclonal Anti-human Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TMN012-0.2MG	0.2mg
TMN012-1MG	1mg

**Description:** Monoclonal Anti-human Neutrophil Gelatinase-associated Lipocalin (NGAL) antibody. This item (TMN012) can be paired with TMN013 for ELISA.

<b>Specificity:</b>	This antibody reacts with human NGAL.
<b>Host Animal:</b>	Balb/c mouse.
<b>Isotype:</b>	IgG1.
<b>Source:</b>	Mouse ascites.
<b>Immunogen:</b>	Recombinant human NGAL.
<b>Format:</b>	Purified, liquid.
<b>Conjugate:</b>	No.
<b>Purification:</b>	DEAE chromatography.
<b>Concentration:</b>	4.25 mg/ml.
<b>Affinity Constant:</b>	Not determined.
<b>Buffer:</b>	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN <sub>3</sub> , pH 7.2.
<b>Pairing Assays:</b>	Pair with TMN013.
<b>Purity:</b>	>95% by HPLC and SDS-PAGE.
<b>Application:</b>	For human NGAL quantitative assays by EIA or CLIA, Lateral Flow, and microfluidics.
<b>Storage:</b>	2-8°C for 1 year. Long-term store at -20°C.

**For Research Use only**