
Specification Sheet

Catalog Number	Unit Size
TBS10101-0.2MG	0.2mg
TBS10101-1MG	1mg

Description:	Mab to clonazepam.
Specificity:	This antibody reacts with clonazepam and 7-clonazepam.
Host Animal:	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse.
Isotype:	IgG1, kappa light chain.
Source:	Mouse ascites.
Immunogen:	Clonazepam-KLH.
Format:	Purified, liquid.
Conjugate:	No.
Purification:	DEAE chromatography.
Concentration:	4.16 mg/ml (OD280 nm, $E^{0.1\%} = 1.40$).
Affinity Constant:	Not determined.
Buffer:	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Not determined.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For clonazepam quantitative assays by EIA.
Storage:	2-8°C for 2 years. Long-term store at -20°C.

For Research Use only