

---

**Specification Sheet**

Catalog Number	Unit Size
TMT061-0.2MG	0.2mg
TMT061-1MG	1mg

<b>Description:</b>	Mab to Tramadol.
<b>Specificity:</b>	Cis-tramadol 100% N-desmethyl-cis-tramadol 70% o-desmethyl-cis-tramadol 0.1%
<b>Host Animal:</b>	Hybridization of P3x63-Ag 9.653 myeloma cells with spleen cells from Balb/c mouse.
<b>Isotype:</b>	IgG2b, kappa light chain.
<b>Source:</b>	Mouse ascites.
<b>Immunogen:</b>	Tramadol-BSA.
<b>Format:</b>	Purified, liquid.
<b>Conjugate:</b>	No.
<b>Purification:</b>	DEAE chromatography.
<b>Concentration:</b>	4.24 mg/ml (OD280 nm, $E^{0.1\%} = 1.40$ ).
<b>Affinity Constant:</b>	Not determined.
<b>Buffer:</b>	0.015 M potassium phosphate buffer, 0.85% NaCl, 0.05% NaN <sub>3</sub> , pH 7.2.
<b>Pairing Assays:</b>	Not determined.
<b>Purity:</b>	>95% by HPLC and SDS-PAGE.
<b>Application:</b>	For tramadol quantitative assays by EIA.
<b>Storage:</b>	2-8°C for 2 years. Long-term store at -20°C.

**For Research Use only**