

Specification Sheet

Catalog Number	Unit Size
TBS10192-100UG	100µg

Description:	Human-Rabbit chimeric antibody is made by fusion of the antigen binding region to rabbit anti-His Tag monoclonal antibody, with the constant domain of human IgG1. The antibody is recognized only by an anti-Human IgG secondary antibody.
Specificity:	This antibody reacts to recombinant proteins containing the 6xHis-Tag or 10xHis-Tag fused to either the amino or carboxy terminus. No cross reactivity with any endogenous protein in mammalian or bacteria cells.
Host Animal:	Rabbit.
Isotype:	IgG.
Immunogen:	Mixture of a peptide with 6xHis-Tag at the N-terminus and another peptide with 6xHis-Tag at the C-terminus.
Format:	Purified, liquid.
Purification:	Protein an affinity purified from an animal origin-free culture supernatant.
Concentration:	1.0 mg/ml.
Unit Size:	100 µg
Buffer:	50% Glycerol/PBS with 0.09% sodium azide.
Application:	For ELISA: Capturer: 0.1 µg/mL – 2 µg/mL; Detector: 0.01 µg/mL – 0.5 µg/mL. For WB: 0.1 µg/mL – 1 µg/mL.
Storage:	Store at -20°C for 1 year.

For Research Use only