

Human Interleukin 6 (IL-6) Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10139H-0.2MG	0.2mg
TBS10139H-1MG	1mg

Description: Human interleukin 6 (IL-6) monoclonal antibody (IL-6 Mab). This item (TBS10139H) is conjugated with HRP. It can be used as a detection antibody to pair with TBS10137 or TBS10138 as a capture antibody for ELISA. The typical data is displayed in Figure 1.

Specificity: This antibody reacts with human IL-6.

Host Animal: Balb/c mouse.

Isotype: IgG2b.

Source: Mouse ascites.

Immunogen: Recombinant human IL-6

Format: Purified, liquid.

Conjugate: HRP

Purification: Chromatography.

Concentration: 1 mg/ml.

Affinity Constant: Not determined.


Buffer: 10 mM PBS, 1% BSA, 50% glycerol, pH 7.2.

Pairing Assays: Paired with TBS10137 or TBS10138.

Purity: >95% by HPLC and SDS-PAGE.

Application: For human IL-6 quantitative assays by EIA or CLIA, Lateral Flow, and microfluidics.

Storage: 2-8°C for 1 year.



For Research Use only