

Phospho-Tau Thr217 (p-tau217) Polyclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10021-0.1MG	0.1mg
TBS10021-0.5MG	0.5mg

Description: Phospho-Tau Thr217 (p-tau217) Polyclonal Antibody (PAB). This antibody (TBS10021) is paired with Mouse Monoclonal Antibody (TBS10022) for ELISA shown in Fig. 1. This antibody can be used for Western Blot (Fig. 2), or IHC in human, rat and mouse brain tissues.

Specificity: This antibody reacts specifically with phosphor- tau Thr217.

Species reactivity: React with: Human, Mouse and Rat.

Host Animal: Rabbit.

Isotype: IgG.

Immunogen: Synthesized phosphor-tau thr217 peptide

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 1.0 mg/ml.

Buffer: 15% Glycerol, 10 mM PBS, 0.05% NaN₃, pH 7.2.

Application: This antibody (TBS10021) is paired with Mouse Monoclonal Antibody (TBS10022) for ELISA, CLIA, Lateral Flow. This antibody can be used for Western Blot, or IHC in human, rat and mouse brain tissues.

Storage: Store at -20°C for 4 years.

Research Use Only

