

Amyloid Beta 1-42 (A β 42) Mouse Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10029-0.1MG	0.1mg
TBS10029-0.5MG	0.5mg

Description: Amyloid Beta 1-42 (A β 42) Mouse Monoclonal Antibody. This antibody as a capture antibody can be paired with TBS10030 as detection.

Specificity: This antibody reacts with Amyloid Beta 1-42 (A β 42).

Host Animal: Balb/c mouse.

Isotype: IgG1.

Immunogen: Amyloid Beta 1-42 (A β 42) peptide.

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 1.0 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN₃, pH7.2.

Pairing Assays: TBS10030

Purity: >95% by HPLC and SDS-PAGE.

Application: For human Amyloid Beta 1-42 (A β 42) by EIA or CLIA and Lateral Flow.

Storage: -20°C for 2 years.

For Research Use only