

## Amyloid Beta 1-40 (A $\beta$ 40) Mouse Monoclonal Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TBS10032-0.1MG	0.1mg
TBS10032-0.5MG	0.5mg

**Description:** Amyloid Beta 1-40 (A $\beta$ 40) Mouse Monoclonal Antibody. This antibody as a detection antibody can be paired with TBS10031 as a capture antibody.

**Specificity:** This antibody reacts with Amyloid Beta 1-40 (A $\beta$ 40).

**Host Animal:** Balb/c mouse.

**Isotype:** IgG1.

**Immunogen:** Amyloid Beta 1-40 (A $\beta$ 40) peptide.

**Format:** Purified, liquid.

**Purification:** DEAE chromatography.

**Concentration:** 1.0 mg/ml.

**Affinity Constant:** Not determined.

**Buffer:** 10 mM PBS, 0.05% NaN<sub>3</sub>, pH7.2.

**Pairing Assays:** TBS10031.

**Purity:** >95% by HPLC and SDS-PAGE.

**Application:** For human Amyloid Beta 1-40 (A $\beta$ 40) by EIA or CLIA and Lateral Flow.

**Storage:** -20°C for 2 years.

**For Research Use only**