

Phospho-Tau Thr231 (p-tau231) Polyclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10035-0.1MG	0.1mg
TBS10035-0.5MG	0.5mg

Description: Phospho-Tau Thr231 (p-tau231) Polyclonal Antibody. This antibody is validated as capture in ELISA.

Specificity: This antibody reacts specifically with phosphor-tau Thr231.

Host Animal: Rabbit.

Isotype: IgG.

Immunogen: Synthesized phosphor-tau thr231 peptide

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 1.0 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN₃, pH7.2.

Purity: >95% by HPLC and SDS-PAGE.

Application: For human phosphor-tau Thr231 quantitative assays by EIA or CLIA and Lateral Flow.

Storage: Store at -20°C for 2 years.

For Research Use only