

Anti-human Beta-2 Microglobulin Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10006-0.2MG	0.2mg
TBS10006-1MG	1mg

Description: Anti-human Beta-2 Microglobulin Mab to B2 microglobulin. This item (TBS10006) can be paired with TBS10005 for ELISA.

Specificity:	This antibody reacts specifically with b2 microglobulin.
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	Human beta-2 microglobulin.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	4.34 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10005 for human beta-2 microglobulin detection.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human beta-2 microglobulin quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	Store at -20°C for 2 years.

For Research Use only