

## Anti-human CA19-9 Monoclonal Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TBS10098-0.2MG	0.2mg
TBS10098-1MG	1mg

**Description:** Anti-human CA19-9 Mab to CA19-9. This item (TBS10098) can be paired with TBS10097 for ELISA.

<b>Specificity:</b>	This antibody reacts specifically with CA19-9.
<b>Host Animal:</b>	Balb/c mouse.
<b>Isotype:</b>	IgG1.
<b>Source:</b>	Mouse ascites.
<b>Immunogen:</b>	CA19-9.
<b>Format:</b>	Purified, liquid.
<b>Purification:</b>	DEAE chromatography.
<b>Concentration:</b>	4.65 mg/ml.
<b>Affinity Constant:</b>	Not determined.
<b>Buffer:</b>	10 mM PBS, 0.05% NaN <sub>3</sub> , pH 7.2.
<b>Pairing Assays:</b>	Pair with TBS10097 for CA19-9 detection.
<b>Purity:</b>	>95% by HPLC and SDS-PAGE.
<b>Application:</b>	For human CA19-9 quantitative assays by EIA or CLIA and Lateral Flow.
<b>Storage:</b>	-20°C for 2 years.

**For Research Use only**