

Anti-human CA15-3 Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10099-0.2MG	0.2mg
TBS10099-1MG	1mg

Description: Anti-human CA15-3 Mab to CA15-3. This item (TBS10099) can be paired with TBS10100 for ELISA.

Specificity: This antibody reacts specifically with CA153.

Host Animal: Balb/c mouse.

Isotype: IgG1.

Source: Mouse ascites.

Immunogen: CA15-3.

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 4.31 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN₃, pH 7.2.

Pairing Assays: Pair with TBS10100 for CA153 detection.

Purity: >95% by HPLC and SDS-PAGE.

Application: For human CA15-3 quantitative assays by EIA or CLIA and Lateral Flow.

Storage: -20°C for 2 years.

For Research Use only