

Anti-human CA-125 Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10104-0.2MG	0.2mg
TBS10104-1MG	1mg

Description: Anti-human CA125 Mab to CA125. This item (TBS10104) can be paired with TBS10103 for ELISA.

Specificity: This antibody reacts specifically with CA125.

Host Animal: Balb/c mouse.

Isotype: IgG1.

Source: Mouse ascites.

Immunogen: CA125.

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 5.02 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN₃, pH 7.2.

Pairing Assays: Pair with TBS10103 for CA125 detection.

Purity: >95% by HPLC and SDS-PAGE.

Application: For human CA125 quantitative assays by EIA or CLIA and Lateral Flow.

Storage: -20°C for 2 years.

For Research Use only