

Anti-human CA72-4 Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10108-0.2MG	0.2mg
TBS10108-1MG	1mg

Description: Anti-human CA72-4 Mab to CA72-4. This item (TBS10108) can be paired with TBS10107 for ELISA.

Specificity:	This antibody reacts specifically with CA72-4.
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	CA72-4.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	4.58 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TBS10107 for CA72-4 detection.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human CA72-4 quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	-20°C for 2 years.

For Research Use only