

Anti-human CA19-9 Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TMC061-0.2MG	0.2mg
TMC061-1MG	1mg

Description: Anti-human CA19-9 Mab to CA19-9. This item (TMC061) can be paired with TMC062 for ELISA.

Specificity: This antibody reacts specifically with CA19-9.

Host Animal: Balb/c mouse.

Isotype: IgG1.

Source: Mouse ascites.

Immunogen: CA19-9.

Format: Purified, liquid.

Purification: DEAE chromatography.

Concentration: 4.93 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 0.05% NaN₃, pH 7.2.

Pairing Assays: Pair with TMC062 for CA19-9 detection.

Purity: >95% by HPLC and SDS-PAGE.

Application: For human CA19-9 quantitative assays by EIA or CLIA and Lateral Flow.

Storage: -20°C for 2 years.

For Research Use only