

Anti-human CA-125 Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TMC092-0.2MG	0.2mg
TMC092-1MG	1mg

Description: Anti-human CA125 Mab to CA125. This item (TMC092) can be paired with TMC091 for ELISA.

Specificity:	This antibody reacts specifically with CA125.
Host Animal:	Balb/c mouse.
Isotype:	IgG1.
Source:	Mouse ascites.
Immunogen:	CA125.
Format:	Purified, liquid.
Purification:	DEAE chromatography.
Concentration:	5.02 mg/ml.
Affinity Constant:	Not determined.
Buffer:	10 mM PBS, 0.05% NaN ₃ , pH 7.2.
Pairing Assays:	Pair with TMC091 for CA125 detection.
Purity:	>95% by HPLC and SDS-PAGE.
Application:	For human CA125 quantitative assays by EIA or CLIA and Lateral Flow.
Storage:	-20°C for 2 years.

For Research Use only