

## Anti-human CA72-4 Monoclonal Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TMC121-0.2MG	0.2mg
TMC121-1MG	1mg

**Description:** Anti-human CA72-4 Mab to CA72-4. This item (TMC121) can be paired with TMC122 for ELISA.

**Specificity:** This antibody reacts specifically with CA72-4.

**Host Animal:** Balb/c mouse.

**Isotype:** IgG1.

**Source:** Mouse ascites.

**Immunogen:** CA72-4.

**Format:** Purified, liquid.

**Purification:** DEAE chromatography.

**Concentration:** 4.96 mg/ml.

**Affinity Constant:** Not determined.

**Buffer:** 10 mM PBS, 0.05% NaN<sub>3</sub>, pH 7.2.

**Pairing Assays:** Pair with TMC122 for CA72-4 detection.

**Purity:** >95% by HPLC and SDS-PAGE.

**Application:** For human CA72-4 quantitative assays by EIA or CLIA and Lateral Flow.

**Storage:** -20°C for 2 years.

**For Research Use only**